
HIGH-FREQUENCY / RADIOSURGICAL DEVICES PRODUCT OVERVIEW

RF Blade* PS-50





- -Monopolar cutting
- -Monopolar coagulation -also via monopolar forceps
- -Monopolar pure coagulation
- -Bipolar coagulation
- -Hand piece or foot switch operation
- -memory of last setting used
- -Switching to the desired output mode by pressing the buttons at the handpieced

RF Blade* PSP-100

Output: maximum 60W/4MHZ

The universal device for all minor and major surgeries making the scalpel superfluous.




- Monopolar cutting
- -Monopolar coagulation -also via monopolar forceps
- -Monopolar pure coagulation
- -Bipolar coagulation
- -Hand piece or foot switch operation
- -memory of last setting used
- -Switching to the desired output mode by pressing the buttons at the handpieced

Output: maximum 120W/2.2MHZ

RF Blade* PSP-100 'P'

the sensational new development in the field of electrosurgery for all minor and major surgeries




- -Monopolar cutting with simultaneous coagulation
- -Monopolar coagulation -also via monopolar forceps
- -Monopolar coagulation with smart switch-off
- -Monopolar spray coagulation
- -Bipolar coagulation
- -Bipolar coagulation with smart switch-on/off
- -Hand piece or foot switch operation
- -4 individual setting can be saved
- -Intelligent tissue resistance detection system
- -Switching to the desired output mode by pressing the buttons at the hand pieced or squeezing the tips of the bipolar forceps together

Output: maximum 120W/2.2MHZ

PARSAM SANAT BEHDAD Pardis Technology Park 20th km of Damavand Road Tehran I.R. Iran









info@psb-co.ir

